



Helm Dental Laboratory

2801 Capital Street
Wylie, TX 75098
972-442-9772

www.HelmDentalLaboratory.com
info@HelmDentalLaboratory.com

ALL information required. Complete all sections. Please keep a copy for your files. We will call you before proceeding if we have any questions.

Restorative Doctor's name

Surgical Doctor's Name (for case reference only)

Address

City State Zip

Phone Doctor's e-mail

Patient name Age Male/Female

Doctor's signature

License number / State issued by Today's date

Final Screw-Retained Prosthesis Packages

Acrylic & Titanium Wraparound

Package includes:

- Rapid Try-In
- Wax Try-In
- Final Screw-Retained Denture
- Model work
*****Analog not included*****
- Bridge screws for final screw-retained denture

Upgrades per arch:

- Additional or replacement bridge screws, each
- Restorative chairside assistance, per visit (e.g impressions, try-in, final seating, etc.)
- **Thermoguard Nightguard**
(highly recommended for protection of prosthesis overnight)

Contact Helm Dental Laboratory for RX with current pricing on this product/package.

All trademarks and company names are the property of their respective owners.
©Helm Dental Laboratory 2-23.

SCREW-RETAINED DENTURE RX FOR TITANIUM & ACRYLIC PROSTHESIS

Case Information

Arch Selection (Please use separate RX for each arch)

- Maxilla
- Mandible

Final Screw-retained Denture type

- Titanium & Acrylic

Please send me the following item:

- Rapid Try-In (Appointment #1)
- Wax Try-In (Appointment #2) Shade: _____
- Additional Wax Try-In (Optional Appointment #2A)
- Final Prosthesis (Delivery Appointment)


See opposite side for summary of appointment steps.
Detailed protocol available upon request.

Notes:

Please list the implant manufacturer, type & diameter (e.g., Astra Tech EV 4.8, Nobel Active RP 4.3, Straumann BL 4.8, etc.) & tooth number. Additional charges for more than 5 implants.

Implant Mfg.	Type & Diameter	Tooth #

Please mark missing teeth/ teeth to be extracted with X on tooth chart below.





Summary of Restorative Protocol Screw-Retained Denture Titanium & Acrylic

This protocol outlines the basics of the restorative appointments for completing the delivery of a Titanium & Acrylic Screw-Retained Denture. ***Detailed protocol available upon request.***

RESTORATIVE PHASE

After surgeon releases patient, Restorative doctor sends Screw-Retained Denture RX for Titanium & Acrylic Prosthesis to lab to order a **Rapid Try-In**. The Rapid Try-In is used as a verification index to take a verified master impression.

Allow 10 working days to receive Rapid Try-In.

Appointment 1 - Rapid Try-In (Master Impression and Verification) - Schedule 60 min.

Rapid Try-In appointment is for verification of implant locations.

Items sent by Lab to complete appointment: *Rapid Try-In*

Tools and parts needed to complete appointment: *Tray adhesive, light body impression material, bite registration material, shade tabs, appropriate torque wrench and screw driver.*

Forms: *Screw-Retained Denture RX for Titanium & Acrylic Prosthesis*

- Complete Rapid Try-In process.
- Complete Screw-Retained Denture RX and select Wax Try-In. Indicate shade for final denture.
- Send RX, Rapid Try-In impression, bite, photos and any scans taken to Lab.

Schedule next appointment after 10 working days in lab.

Appointment 2 - Wax Try-In - Schedule 60 min.

Try-in appointment is to verify esthetics and function for final restoration, including final shade.

Items sent by Lab to complete appointment: *Wax Try-In*

Tools & parts needed to complete appointment: *Light body impression material, shade tabs, appropriate torque wrench and screw driver.*

Forms: *Screw-Retained Denture RX, Try-In Checklist, Process & Finish Form.*

- **If there are significant changes** in tooth placement or fit, another Wax Try-In is recommended. This will add an extra appointment to the process. Send RX, Wax Try-In, and Try-In Checklist to Lab.
- **If there are not significant changes**, send RX, Wax Try-In, Try-in Checklist, and signed Process & Finish Form to Lab. Any changes requested on the finished prosthesis after approval of the Try-In will incur additional charges.

Schedule next appointment after 15 working days in lab.

Appointment 3 - Final Prosthesis Delivery - Schedule 30-45 min

Items sent by Lab to complete appointment: *Final Titanium & Acrylic Screw-Retained Prosthesis*

Tools & parts needed to complete appointment: *Light body impression material, Teflon tape, composite, appropriate torque wrench and screw driver.*

Forms: *None*

- Seat **Final Prosthesis**.